SprintRay EU Temporary Crown & Teeth Instruction for use

1. Indication

SprintRay EU Temporary Crown & Teeth is a resin for the generative production of long term temporary dental restorations based on image projection systems.

2. Contraindication

SprintRay EU Temporary Crown & Teeth restorations are contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... for bridges with more than one pontic.
3. ... for every application that is not part of the indication (see above).

3. Requirements

Software:

- RayWare (CAD)

Hardware:

- Pro Desktop 3D (3D Printing)
- Pro Wash / Dry (Cleaning)
- ProCure (Post Curing)

SprintRay, Inc., 3577 North Figueroa St., Los Angeles, CA 90065

4. Material

SprintRay EU Temporary Crown & Teeth consists of functional (meth)acrylic resins and inorganic fillers with particle sizes from 0.4 to 3 microns.

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occlusal 1.5 mm (central fissure)
circumferential 1 mm

6. Material parameters

Radiation penetration depth controllable by exposure time

≥ 100 µm (4 mils)
≥ 125 µm (5 mils)

7. Manufacturing process (fig. 1-10)

1. Prepare data (RayWare).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (ProWash / Dry):
- blow surplus resin out of cavities with pressured air before washing
- wash 5 min. / 5 min.
8. Dry parts (ProWash / Dry):
- 3 min.
- additional: dry with pressured air as required
9. Post curing (ProCure): 5 min. / 60°C.
10. Finish parts.

8. Individualization and cementation

The temporary restorations can be individualized by light curing composite materials. The reduction of the restoration should not exceed 0.3 mm (incisal and vestibular area). Please follow the instructions of the light curing composite manufacturers. For cementation common provisional dental cements (free of Eugenol) and adhesives are suitable. The inner areas of the crown should be conditioned by sandblasting or other commonly used techniques. Please follow the instructions and indications of the corresponding manufacturers.

9. Notice


Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.

Caution: Polymerized resins are chemically resistant - avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

10. Hazard statements

Causes skin irritation. May cause an allergic skin reaction.
Causes serious eye damage. May cause respiratory irritation.
Harmful to aquatic life with long lasting effects.

UMDNS 16-697

Physical properties*/
Physikalische Eigenschaften*/
Propriétés physiques*/
Propiedades físicas*/
Proprietà fisiche*:

SprintRay EU Temporary Crown & Teeth

■ Flexural strength/
Biegefestigkeit/
Résistance à la flexion/
Fuerza flexible/
Resistenza alla flessione:
MPa (ISO 4049):
> 100

■ Modulus of elasticity/
Elastizitätsmodul/
Module d'élasticité/
Modulo di elasticità/
Módulo de elasticidad:
MPa (ISO 4049):
> 2000

■ Water absorption/
Wasseraufnahme/
Absorción de'eau/
Assorbimento dell'acqua:
complies with/ erfüllt/
est conforme/ conforme a/
cumple con ISO 10477

■ Solubility/
Löslichkeit/
Solubilité/
Solubilidad/
Solubilità:
complies with/
est conforme/
cumple con/
conforme a
ISO 10477

■ Inorganic filler/
Anorganischer Füllstoff/
Charge inorganique/
Relleno inorgánico/
Riempitivo inorganico:
40.0 m-%

Ordering information/
Bestellinformationen/
Informations de commande/
Información sobre pedidos/
Informazioni sull'ordine:

SprintRay EU Temporary Crown & Teeth

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