# SprintRay EU Castable Instruction for use

## 1. Indication

**SprintRay EU Castable** is a resin for the generative production of dental castforms based on image projection systems.

## 2. Contraindication

**SprintRay EU Castable** is contraindicated for every application that is not part of the indication (see above).

## 3. Requirements

### Software:

- RayWare (CAD)

### Hardware:

- Pro Desktop 3D (3D Printing)  
- Pro Wash / Dry (Cleaning)  
- ProCure (Post Curing)

SprintRay, Inc., 3577 North Figueroa St., Los Angeles, CA 90065

## 4. Material

**SprintRay EU Castable** consists of functional (meth)acrylic resins.

## 5. Geometric presetting

According to the software design specifications.

## 6. Material parameters

Radiation penetration depth controllable by exposure time

25 µm  
50 µm  
100 µm

## 7. Manufacturing process (fig. 1-10)

1. Prepare data (RayWare).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (ProWash / Dry): 5 min. / 5 min.
8. Dry parts (ProWash / Dry): 3 min.
9. Post curing (ProCure): 5 min. / 60 °C.
10. Finish parts.

## 8. Finishing processes

Build up of material or repair possible with **LCR-1** gel.

## 9. Notice


Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/ printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.

**Caution:** Polymerized resins are chemically resistant - avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

## 10. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. May cause long lasting harmful effects to aquatic life.

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