# SprintRay EU Splint Flex Instruction for use

## 1. Indication

**SprintRay EU Splint Flex** is intended for the production of orthodontic and dental objects (mouthguards, nightguards, snoring appliances and splints).

## 2. Contraindication

**SprintRay EU Splint Flex** is contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... for denture bases.
3. ... for every application that is not part of the indication (see above).

## 3. Requirements

### Software:

- RayWare (CAD)

### Hardware:

- Pro Desktop 3D (3D Printing)
- Pro Wash / Dry (Cleaning)
- ProCure (Post Curing)

SprintRay, Inc., 3577 North Figueroa St., Los Angeles, CA 90065

## 4. Material

**SprintRay EU Splint Flex** consists of functional (meth)acrylic resins.

## 5. Geometric presetting

Minimum wall thickness: 1.5 mm

## 6. Material parameters

Radiation penetration depth controllable by exposure time  
50 µm  
100 µm

## 7. Manufacturing process (fig. 1-10)

1. Prepare data (RayWare).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (ProWash / Dry): 5 min. / 5 min.
8. Dry parts (ProWash / Dry): 3 min.
9. Post curing (ProCure): 20 min. / 60 °C.
10. Finish parts.

## 8. Finishing processes

Polishing

## 9. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/ printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.


According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

**Caution:** Polymerized resins are chemically resistant - avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

## 10. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Harmful to aquatic life with long lasting effects.

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
**Software:**  
- RayWare (CAD)

### Hardware:

- Pro Desktop 3D (3D Printing)  
- Pro Wash / Dry (Cleaning)  
- ProCure (Post Curing)

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