




Study Model White

Rapidly 3D Print Study Models for Patient Presentation

SprintRay Study Model White is a special material formulation that provides outstanding color, clarity, and smoothness, making it perfect for patient presentation. With almost no pigment separation, you can print Study Model White with confidence, knowing that your models will have fantastic color distribution and results that are reliable and accurate.

SprintRay Study Model White features improved flexural strength and modulus for a sturdier model that is less likely to break or deform. With better platform adhesion and reduced odor, SprintRay Study Model White can help take your practice to the next level, delighting patients with striking visual presentations that help improve patient education and case acceptance.



SRI-0202008

SKU Number

KEY BENEFITS

- ✓ Improved flexural modulus and strength
- ✓ Up to 10% faster print speeds
- ✓ Print four full arch study models in just 20 minutes*
- ✓ Deliver outstanding patient experience by bringing results to life
- ✓ Rich color and smooth surface mimic orthodontic models

* When printed on SprintRay Pro with Ludicrous Speed.

FLEXURAL STRENGTH

Higher flexural strength ensures that the model doesn't break during normal use.

Study Model White 91.5 MPa

Die and Model Tan 66.7 MPa

FLEXURAL MODULUS

Higher flexural modulus prevents the model from distorting, which can impact accuracy.

Study Model White 2650 MPa

Die and Model Tan 1700 MPa