



THE CHAIRSIDE 3D PRINTED

**CROWN**

---

PORTFOLIO

# SPRINTRAY CROWN 3D PRINTING. A COMPLETE CHAIRSIDE SOLUTION.

SprintRay Crown is a ceramic filled hybrid material for 3D printing single crowns, inlays, onlay, and veneers. This material is easy to finish and polish using traditional methods, and features excellent aesthetics to blend with existing teeth.

## A New Class Of Resin, Unlike Any Other.

High strength and resistance to wear make SprintRay Crown a go-to material for restorative 3D printing. And because SprintRay Crown fulfills all the requirements for a Class II medical device, you can place with confidence.


Crown is available in seven shades:  
A1, A2, A3, B1, B3, C2, and D3.



## Key Features

- Biocompatible solution that fulfills all the requirements for a Class II medical device
- Versatile: use as a single crown on natural tooth or abutment, as an inlay, onlay or veneer
- Very low water absorption reduces the tendency to age while a smooth surface means low plaque accumulation

# STREAMLINE EVERYTHING. PRINT DEFINITELY.




## AI-Powered Crown Design

Let powerful artificial intelligence be your expert guide – AI automates the design process and makes it completely hands-free. Posterior crown designs are turned around within minutes, allowing for single-visit restorative dentistry.



Posterior full contour crown.  
Printed at 100  $\mu$ m layer thickness.



## Advanced 3D Printing

The Pro S Crown Kit consists of a new compact build platform and resin tank engineered for rapid production of up to 6 single-unit crowns from scan to delivery in less than 45 minutes.



## Post Processing Made Easy

Featuring ProCure 2, the fastest, most advanced post-curing system on the planet. Cure printed crowns in just a few minutes.



T H E C O M P L E T E C H A I R S I D E S O L U T I O N

# SPINTRAY CROWN 3D PRINTING. A COMPLETE CHAIRSIDE SOLUTION.



## Proven Materials

Experience a new class of ceramic dominant, Class II medical device resin, designed for 3D printing definitive restorations.

## Delivery in Less than 45 Minutes

Increase case acceptance by offering same-day scan-to-delivery thanks to AI design and the Pro S Crown Kit.

## Cost-efficient

A complete solution set to transform same-day, chairside delivery for a price you won't believe.




SPINTRAY. THE ONLY COMPLETE 3D PRINTING SOLUTION IN DENTISTRY



SprintRay Europe GmbH  
Brunnenweg 11  
64331 Weiterstadt  
Germany

[www.sprintray.com](http://www.sprintray.com)  
[info.eu@sprintray.com](mailto:info.eu@sprintray.com)



Learn More

